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MEISINGER Digital Surgical guides

EN

MEISINGER surgical guides are individual custom-made products produced by the dental technician according to the intended prosthetic treatment result. They reflect the position of the implants on the alveolar process in metric relation to the subsequent prosthetic restoration. Thus, they support the implementation of implant planning and enable the exact placement of the implant at the planned position in the patient's jaw. They can be used independently of the implant system.

Indications for use

MEISINGER surgical guides are indicated for positioning of implants, as well as for the guidance of surgical instruments for the implant site preparation and the corresponding implants during insertion.

Material description

NextDent™ SG (Surgical Guide, manufacturer: Vertex-Dental B.V., Netherlands) is a biocompatible Class I material, developed for the printing of Surgical Guides for implant surgery use. Because of the high precision of this material, it is easy to insert drill sleeves, directly after printing. Enabling even greater precision during surgery. The properties of NextDent™ SG will not change by using disinfectants. In addition, the material can also be sterilized using standard autoclave protocols. The use of an autoclave does not affect the dimensional stability, therefore NextDent™ SG can be used in every operating theatre.

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Plastic and packaging waste

- The product NextDent™ SG in its polymerized form is not environmentally harmful.

Cleaning and sterilization instructions

- NextDent™ 3D-printing material should be cleaned with nonchemical products. If disinfecting before intended use is required, an ethanol solution can be used.
- NextDent™ SG can be sterilized by making use of an autoclave or using gamma-ray sterilization.
- Do not use a dishwasher or thermal disinfectant.

NextDent™ recommends the following method of sterilization:

- Place the surgical guide in a standard sterilization pouch
- Apply sterilization by placing the pouch containing the guide in an autoclave for 15 minutes at 121°C or 3 minutes at 138°C.
 - Make sure no mechanical forces are applied to the guide during sterilization.
- Let the guide cool down to room temperature before using the guide.
 - Make sure no mechanical forces are applied to the guide during cooling down.

Contraindications

In principle, general medical as well as local, absolute and relative contraindications must be observed for MEISINGER surgical guides.

- Known allergies or hypersensitivities to the chemical components of NextDent™ SG
- Insufficient tooth status for stable support of the surgical guide

General instructions

The MEISINGER Digital surgical guide must only be used by dentists, doctors and / or respective experts who, due to their training and experience, are intensely familiar with the use of this product and who have the corresponding expertise. It is the sole responsibility of the doctor in charge who, depending on the respective situation (indication) decides on the actual use of this product for each individual case. The product is intended for single use only.

The surgical guide must be processed before first use on the patient. Inappropriate processing can result in the patient being infected with harmful bacteria.

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Non-sterile

Follow the instruction for use

Under U.S. federal law, this device may only be sold to or on behalf of trained medical professionals.

Do not use if the packaging is damaged

**Hager & Meisinger GmbH**

Hansemannstr. 10
41468 Neuss • Germany
Tel. +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de

**Hager & Meisinger GmbH**

Hansemannstr. 10
41468 Neuss • Germany
Tel. +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de