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Application and safety instructions for the MEISINGER

## Periimplantitis Kit

### Tungsten Carbide Finishers for Periimplantitis Treatment

Art.-No.: 2575

The Periimplantitis Kit contains eight perfectly matched tungsten carbide finishing burs in egg and flame shapes for the ideal peri-implantitis treatment. Both shapes come in two sizes and each in standard and ultra-fine toothing. The shapes, sizes and toothings of the working parts are perfectly suited for the intraoral treatment of titanium and for different shapes of implant neck and shoulder. All finishers have an extra-long FG shank (total length 32 mm) which allows particularly the treatment of difficult-to-access deep areas. For the peri-implantitis treatment, the finishers are used clockwise and are guided around the exposed implant counterclockwise so the implant surface is smoothed.

## Indications

The Bone Management® Periimplantitis Kit is indicated for the intraoral processing and smoothing of exposed titanium implant surfaces for the surgical treatment of peri-implantitis.

## Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

#### Preparation

The initial situation is an exposed titanium implant surface after bone resorption due to existing peri-implantitis. As preparation for the optimal peri-implantitis treatment, the prosthetics have to be removed from the implant. The implant is provided with a native cover screw, so that the interface is protected against contamination.

#### Processing of the titanium implant surface

Prior to the processing of the titanium surface, the contaminated tissue around the implant is removed with the aid of an appropriate manual instrument (e.g. Curette Art.-No. HI317). Then, the surface structure of the implant is removed with the fine tungsten carbide finishers before the titanium is smoothed with the ultra-fine finishers. During the processing **the clockwise rotating instruments are led around the exposed surface counterclockwise** and with low contact pressure, without damaging the implant shoulder. This kind of use allows for the best possible control of the used instruments. Size and shape of the instruments depend on the existing patient situation and the implant shape.

When using the ultra-fine finishing bur a very smooth surface structure is achieved. The use of elastic polishers can and should be avoided, also to avoid silicone residues.

#### Post-processing

Prior to the tension-free wound closure all metal residues have to be removed by sufficient rinse and/or the application of appropriate prophylaxis brushes (e.g. Art.-No. 115, 116).

#### Note

During the whole processing time, attention has to be paid to sufficient cooling with sterile cooling liquid. Furthermore, the implant shoulder may not be damaged, so that the implant-abutment-interface remains intact also after the treatment. Excessive pressure during the use of the instruments has to be avoided. Also the use of the instruments with leverage effect is not allowed, due to an increased risk of breakage.

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**Contraindications**

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.


**Adverse reaction**

By smoothing the implant surface, contamination with metal shavings can not be completely avoided. This may possibly promote wound healing disorders.

**NOTE:** This procedure is a surgical therapy. Therefore, all alternative, non-surgical treatment options should be considered in advance and the best possible procedure, according to available data, should be selected. If the progression of the inflammation has already progressed too far or should not be stopped in the treatment, it should be assessed whether an explantation of the implant should be performed in order to avoid further bone resorption.

**General instructions**

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

**Hager & Meisinger GmbH**  
Hansemannstr. 10  
41468 Neuss • Germany  
Tel. +49 2131 2012-0  
Fax: +49 2131 2012-222  
Email: info@meisinger.de  
Internet: www.meisinger.de

Distributed by: **Meisinger USA, L.L.C.**  
10150 E. Easter Avenue  
Centennial, Colorado 80112 • USA  
Phone: +1 (303) 268-5400  
Fax: +1 (303) 268-5407  
Email: info@meisingerusa.com  
Internet: www.meisingerusa.com

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