


*my*plant  
two

BACK TO  
THE ROOTS,  
SHAPING THE  
FUTURE

Made in Germany. Made by MEISINGER.

PRODUCT  
CATALOGUE





## THE TWO-PIECE IMPLANT SYSTEM

The concept of a two-piece titanium implant with special expansion thread and self-locking cone as abutment connection was already developed by Prof. Dr. med. dent. Nentwig and Dr. Ing. Moser 30 years ago and brought to clinical maturity. High primary stability, minimum construction height, micro-movement-free and a bacteria-proof implant-abutment connection as well as deep platform switching soon proved to a superior combination in terms of achievable bone and soft tissue stability, and thus guarantors of exceptional long-term success.

With the **myplant two** implant system, Prof. Nentwig and Dr. Moser have further advanced and optimized this concept, which has been documented over decades, and adapted it to the criteria of modern, future-oriented implant therapy.

We, as myplant GmbH, have made it our mission to provide patients and users of **myplant two** with an implant system that provides the best possible preconditions to ensure long-term implant success with sustainable hard and soft tissue stability.

With more than 30 years of dental implant manufacturing and development competence, as well as many years of experience in the commercialisation of implants, myplant GmbH represents an alliance that gives the user the reassuring feeling that **myplant two** is a long-term functional, mechanically stable, as well as tissue-compatible and sophisticated aesthetic solution.

**myplant**  
**two**

Made in Germany. Made by



# CONTENT

## SYSTEM

The myplant two system	2-7
The prosthetic concept	8-11

---

## PRODUCTS


Implants	12
Prosthetic components	13-24
Surgical instruments	25-30
Accessories	31-33
Surgical kits	34-37
Prosthetic kits	38-41

---

## FURTHER INFORMATION

Quick overview of surgical procedures	42-43
Recommended torques for myplant two prosthetics	44
Material specifications	45
Important notes on our ordering services	46
myplant two implant guarantees	47-50
Index	51-56

---



# SYSTEM

## THE MYPLANT TWO SYSTEM

The myplant two system is based on the proven principles of over 30 years in terms of progressive thread design and a bacteria-proof, self-locking, rotation-stable conical connection.



Deep **platform switching** provides an increased area for bone deposition (on the implant shoulder)

Can be inserted **subcrestally** by up to 2 mm if bone supply is adequate

**Self-locking, absolutely rotation-stable internal conical connection** provides a virtually bacteria-proof seal

**Apical bevel** for simplified insertion of the implant

Free positioning of angled abutments through **non-indexed conical connection**

**Rounded implant tip** for gentle sinus floor elevation






# THE IMPLANTS

The myplant two implant system is a two-piece system based on state-of-the-art technology. The implants are manufactured from Grade 4 pure titanium according to ISO 5832-2. This material stands for the highest level of biocompatibility and outstanding mechanical properties, thus offering optimal conditions for secure osseointegration. Grade 4 pure titanium demonstrates a perfect combination of ductility and strength. It offers excellent corrosion resistance and does not contain any toxic constituents. Owing to its outstanding properties, Grade 4 pure titanium has been successfully employed as a material for implantology and general medical technology for over 25 years.

myplant two implants are available in three diameters and five different lengths. Due to the practice-oriented graduation of implant sizes the system is suitable for all indications in dental implantology, even in difficult bone conditions.

The letter and color coding system allows fast and safe identification of the various implant diameters. Corresponding instruments for the implant bed preparation are marked with the same color code.

## Color coding

	Red	Implant diameter 3.5 mm
	Orange	Implant diameter 4.0 mm
	Yellow	Implant diameter 4.5 mm
















The implant name includes a capital letter which, same as the colour, identifies the implant diameter. The following numbers define the length of the implant in millimeters.



L [mm]					
	6.6	8.0	9.5	11.0	14.0
Ø [mm]					
3.5	A 6.6	A 8	A 9.5	A 11	A 14
4.0	M 6.6	M 8	M 9.5	M 11	M 14
4.5	B 6.6	B 8	B 9.5	B 11	B 14



## INDICATIONS AT A GLANCE

IMPLANTS			PROSTHETICS		
myplant two	Length	Anterior   Canine   Premolar   Molar	Single Tooth	Bridge	Telescope Crowns   Ball Anchors   LOCATORS®
Ø 3.5 mm	 6.6 mm	✓*	X	✓	✓*
	 8.0 mm	✓	✓	✓	✓
	 9.5 mm	✓	✓	✓	✓
	 11.0 mm	✓	✓	✓	✓
	 14.0 mm	✓	✓	✓	✓
Ø 4.0 mm	 6.6 mm	✓*	X	✓	✓*
	 8.0 mm	✓	✓	✓	✓
	 9.5 mm	✓	✓	✓	✓
	 11.0 mm	✓	✓	✓	✓
	 14.0 mm	✓	✓	✓	✓
Ø 4.5 mm	 6.6 mm	✓*	X	✓	✓*
	 8.0 mm	✓	✓	✓	✓
	 9.5 mm	✓	✓	✓	✓
	 11.0 mm	✓	✓	✓	✓
	 14.0 mm	✓	✓	✓	✓



Special indications for implants with 6.6 mm length:

Due to the reduced surface for anchorage in the bone, the implants with a length of 6.6 mm are only recommended for the following indications until appropriate clinical studies are available:

- In edentulous jaws: as auxiliary implant/supporting implant for implant-supported bar constructions or splinted bridges.
- Partially edentulous jaws: as auxiliary implant/supporting implant in connection with longer implants splinted to a suprastructure.

**Important:** Attention should be paid to the load distribution of the restoration.

## THE MACRO DESIGN


### The thread design

The progressive thread design in combination with the three-stage preparation technique lead to very high primary stability, even in cases of average bone qualities. The thread depth increases in apical direction and shows an arch-shaped curve of the thread shanks, thus achieving favorable biomechanical load distribution into the bone. The apically enlarging thread depth allows good anchorage in differing bone qualities as well as promoting a bone stimulating load distribution during mastication. Vertical and lateral forces are primarily deflected to the elastic cancellous bone, whereas the cortical bone is relieved, which is essential for the long-term maintenance of the marginal bone level and the resulting esthetics. Stress concentration in the emergence area of the implant, as proven for implants with consistent threads, is thus avoided. The preparation technique and the special macro design are matched perfectly to the natural bone structure and result in high primary stability with maximum bone to implant contact, even when bone quality is compromised.



### The implant surface


The enossal surface of the implants is blasted with corundum and thus creates a macro-roughness on the titanium surface. This is followed by acid etching adding a micro-roughness to the implant surface. The resulting maximum increase in surface area promotes the ongrowth of bone tissue, leading to a stable implant-bone connection and supporting the natural healing process. In contrast to many other systems, this surface treatment is also performed on the implant shoulder. Therefore a subcrestal implant insertion is possible which reduces stress during the healing period enabling a strong osseointegration and consequently offering superior support of peri-implant tissues.





## SUBCRESTAL PLACEMENT

In contrast to most implant systems, the myplant two implant was developed specifically for subcrestal placement and can be inserted up to 2 mm below bone level if bone supply is adequate and taking into account the surrounding structures. This allows the crestal bone to grow over the implant shoulder up to the abutment diameter, which benefits increased implant stability and biological support of the peri-implant soft tissue. Various cover screws are available to avoid complete overgrowth of the implant with bone during submerged healing. Each implant packaging includes a sterile cover screw exceeding the implant by 1 mm.



© Prof. Dr. med. dent. Georg-Hubertus Nentwig



**Note:**

Subcrestal placement must be taken into consideration when selecting the length of the implant.

If epicrestal implant positioning is desired or indicated, a flush cover screw is available separately.

## Stable soft tissue support

Deep platform switching and the resulting wide implant shoulder allow more interproximal space at the abutment level than non-conical connections. In combination with the bony deposition on the implant shoulder, this is decisive for establishing a stable and healthy soft tissue cuff and consequently for esthetics. Due to the enlarged interproximal space esthetically pleasing results can also be achieved with tightly placed implants.



## THE PROSTHETIC CONCEPT

The key to successful prosthetics is a firm and tight tapered conical connection. The 360° rotation option of the prosthetic components ensures optimum positioning of angled abutments without making any compromise.

A major advantage of the interface is the fact that all implants have the same internal geometry, allowing each prosthetic component to fit into each implant. The choice of implant is made exclusively on the basis of the available bone and is not restricted by the prosthetics. This also keeps storage space and costs as low as possible. All indications, ranging from single crowns via bridges to partial dentures and implant-supported full dentures, can be restored with the various abutments.

**Note:**

The abutment screw is not removable from the (abutments.)



The different abutment series allow friction-based, screw-retained, bonded or cemented fixation to the corresponding abutments.

To enable cement-free mounting, both the straight as well as the angled abutments are equipped with occlusal screw channels. This not only allows improved esthetics due to deeper crown margins but also avoids the risk of cementitis.



A further alternative in conjunction with standard abutments is the friction-based retention of the superstructure via precisely fitting taper caps, which can be incorporated into the prosthesis without difficulties. Whether polymerized into acrylic or bonded into a tertiary structure, conical caps offer full flexibility.

The shoulder abutments offer the option of cementing or screw-retaining the prosthetic restoration. This offers the clinician the freedom to choose the attachment mode according to need and indication.

When fabricating individual all-ceramic restorations on the titanium base, the fabricated stump is bonded extraorally to the titanium base and the final restoration is fixated in the patient's mouth.

## ABUTMENT REMOVAL


The self-locking cone enables an absolutely rotation-stable implant-abutment-connection which provides a virtually bacteria-proof seal. This connection can be disengaged again with the aid of the abutment remover.

First, the straining screw of the abutment is completely loosened with the screwdriver. The abutment remover is then screwed into the occlusal screw seat. Once resistance is felt, the abutment remover has made contact with the head of the abutment. If it is now rotated further with care, the abutment remover lifts the abutment out of the implant enabling it to be removed easily.

**Note:**

If the abutment screw is not loosened completely, the abutment remover could be blocked.

Should the abutment not come out easily the abutment remover should be disconnected and the screw driver used to fully loosen the straining screw and remove the abutment from the implant.




**Important:** To avoid jamming, it is essential to ensure that the abutment remover is screwed into the screw channel of the occlusal screw seat in the case of angled abutments.

## HIGHLY RESILIENT ABUTMENT-IMPLANT CONNECTION

The special design of the myplant two system results in an outstanding mechanical load capacity with high fatigue strength. The design of the inner cone results in a highly resilient and virtually bacteria-proof connection between abutment and implant.

### Significant stress reduction in the abutment

In combination with different implants



Load according to ISO14801 / 250 N

### Standardized conical connection


The non-indexed conical connection allows free positioning as well as simple parallel alignment of angled abutments.

Despite different implant diameters, all implants have a uniform prosthetic interface. As a result, the choice of implants can be made exclusively depending on the bone supply, while in the selection of the abutments, the focus can be completely placed on the prosthetic requirements. In addition, this results in a reduced prosthetic assortment and leads to lower storage requirements.

## OPTIONS FOR PROSTHETIC RESTORATIONS

All restorations from single crowns via bridges to partial and full dentures can be tackled with the myplant two abutment portfolio. A major advantage is the unified abutment interface of all myplant two implants allowing each prosthetic component to fit into each implant. The different abutment types allow for a friction-retained, screwed, bonded or cemented attachment of the restorations.






Abutmenttype	Single Tooth	Bridge	Partial prosthesis	Total Prosthesis
Standard Abutment	✓	✓	✓	✓
Shoulder Abutment	✓	✓	-	-
Titanium Base	✓	✓	-	-
Multi Unit Abutment	-	✓	✓	✓
LOCATOR®	-	-	✓	✓
Ball Anchor	-	-	✓	✓








# PRODUCTS

The drilling speeds and torques of the instruments and system components are summarised in the overviews on pages 43 and 44.






## IMPLANTS

					
	A6.6	A8	A9.5	A11	A14
Ø (mm)	3.5	3.5	3.5	3.5	3.5
Length (mm)	6.6	8.0	9.5	11.0	14.0
Order No.	A3566	A3580	A3595	A3511	A3514


  

					
	M6.6	M8	M9.5	M11	M14
Ø (mm)	4.0	4.0	4.0	4.0	4.0
Length (mm)	6.6	8.0	9.5	11.0	14.0
Order No.	M4066	M4080	M4095	M4011	M4014

					
	B6.6	B8	B9.5	B11	B14
Ø (mm)	4.5	4.5	4.5	4.5	4.5
Length (mm)	6.6	8.0	9.5	11.0	14.0
Order No.	B4566	B4580	B4595	B4511	B4514

A sterile 1.0 mm high cover screw (Order No. PVS01) is supplied with all implants.



# PROSTHETIC COMPONENTS


## Cover Screws

			
Description	Cover Screw 0 mm	Cover Screw 1 mm	Cover Screw 2 mm
Order No.	PVS00	PVS01	PVS02

## Healing Abutments

			
Description	Healing Abutment Ø 4.2 mm GH 1.5 mm	Healing Abutment Ø 4.2 mm GH 3.0 mm	Healing Abutment Ø 4.2 mm GH 4.5 mm
Order No.	PGF15	PGF30	PGF45

**Note:** The abutment interface is identical for all myplant two implants. As a result, the listed cover screws and gingiva formers can be used for all implants.





## Transfer Post



Description	Transfer Post, S open tray PUP01	Transfer Post, L open tray PUP02
Order No.		

The transfer posts are tightened by hand.

## Repositioning Post



Description	Repositioning Post closed tray PRP01	Repositioning Post, narrow closed tray PRP02
Order No.		

The repositioning post is tightened by hand or with the aid of a screwdriver.







## Implant Analog



Description	Implant Analog
Order No.	PLI01

# STANDARD ABUTMENTS

## Abutments - straight

			
Description	Abutment 0° GH 1.5 mm, H 4.0 mm PAB01	Abutment 0° GH 3.0 mm, H 4.0 mm PAB03	Abutment 0° GH 4.5 mm, H 4.0 mm PAB05
Order No.			
			
Description	Abutment 0° GH 1.5 mm, H 6.0 mm PAB02	Abutment 0° GH 3.0 mm, H 6.0 mm PAB04	Abutment 0° GH 4.5 mm, H 6.0 mm PAB06
Order No.			

### Note:

To enable cement-free mounting, both the straight as well as the angled abutments are equipped with occlusal screw seats.





### Note:

With sufficient primary stability, a temporary immediate restoration can be achieved on the standard abutments with the aid of the healing caps. This requires a primary stability of >30 Ncm!



### Note:

The abutment screw is not removable from the abutments.

## Abutments - angled

		
Description	Abutment 15° GH 1.5 mm, H 4.0 mm PAB51	Abutment 15° GH 3.0 mm, H 4.0 mm PAB53
Order No.		
		
Description	Abutment 15° GH 1.5 mm, H 6.0 mm PAB52	Abutment 15° GH 3.0 mm, H 6.0 mm PAB54
Order No.		



## Transfer Caps

		
Description	Transfer Cap H 4.0 mm	Transfer Cap H 6.0 mm
Order No.	PAK01	PAK02

## Healing Caps

### Note:





The healing caps are used for the fabrication of temporary immediate restorations. Care must be taken to ensure that the temporary restorations are not exposed to occlusal stress.

		
Description	Healing Cap H 4.0 mm	Healing Cap H 6.0 mm
Order No.	PHK01	PHK02



## Laboratory Analog one-piece

				
Description	Laboratory Analog one-piece 0° H 4.0 mm	Laboratory Analog one-piece 0° H 6.0 mm	Laboratory Analog one-piece 15° H 4.0 mm	Laboratory Analog one-piece 15° H 6.0 mm
Order No.	PLA04	PLA06	PLA54	PLA56

## Wax-Up Copings for straight and angled Abutments

				
Description	Wax-Up Coping H 4.0 mm PMK01	Wax-Up Coping H 6.0 mm PMK02	Wax-Up Coping Rotation Secured H 4.0 mm PMK03	Wax-Up Coping Rotation Secured H 6.0 mm PMK04
Order No.				

## Tapered Caps for straight and angled Abutments

				
Description	Taper Cap without Retention H 4.0 mm PKK04	Taper Cap without Retention H 6.0 mm PKK06	Taper Cap with Retention H 4.0 mm PKK14	Taper Cap with Retention H 6.0 mm PKK16
Order No.				

## Occlusal Screw

	
Description	Occlusal Screw
Order No.	POS01

## Parallel Post

	
Description	Parallel Post for angled Abutments
Order No.	ZPS01

### Note:





Connect the parallel posts to the occlusal screw seats of angled abutments to facilitate a much easier parallel alignment between multiple abutments.

## TRIAL ABUTMENTS



Trial Abutment Kit (Order No. PAP00)  
contains the following articles:



### Trial Abutments - straight

				
Description	Trial Abutment 0° GH 1.5 mm, H 4.0 mm	Trial Abutment 0° GH 1.5 mm, H 6.0 mm	Trial Abutment 0° GH 3.0 mm, H 4.0 mm	Trial Abutment 0° GH 3.0 mm, H 6.0 mm
Quantity / Kit	4	4	4	4
Order No.	PAP01	PAP02	PAP03	PAP04





		
Description	Trial Abutment 0° GH 4.5 mm, H 4.0 mm	Trial Abutment 0° GH 4.5 mm, H 6.0 mm
Quantity / Kit	4	4
Order No.	PAP05	PAP06





### Trial Abutments - angled

				
Description	Trial Abutment 15° GH 1.5 mm, H 4.0 mm	Trial Abutment 15° GH 1.5 mm, H 6.0 mm	Trial Abutment 15° GH 3.0 mm, H 4.0 mm	Trial Abutment 15° GH 3.0 mm, H 6.0 mm
Quantity / Kit	4	4	4	4
Order No.	PAP51	PAP52	PAP53	PAP54





# SHOULDER ABUTMENTS





## Abutments - straight

				
Description	Shoulder Abutment 0° Ø 5.7 mm GH 1.5 mm, H 6.0 mm PSA01	Shoulder Abutment 0° Ø 5.7 mm GH 3.0 mm, H 6.0 mm PSA02	Shoulder Abutment 0° Ø 5.7 mm GH 1.5 mm, H 4.0 mm PSA05	Shoulder Abutment 0° Ø 5.7 mm GH 3.0 mm, H 4.0 mm PSA06
Order No.				

				
Description	Shoulder Abutment 0° Ø 4.3 mm GH 1.5 mm, H 6.0 mm PSA11	Shoulder Abutment 0° Ø 4.3 mm GH 3.0 mm, H 6.0 mm PSA12	Shoulder Abutment 0° Ø 4.3 mm GH 1.5 mm, H 4.0 mm PSA15	Shoulder Abutment 0° Ø 4.3 mm GH 3.0 mm, H 4.0 mm PSA16
Order No.				

## Abutments - angled


				
Description	Shoulder Abutment 15° Ø 5.7 mm GH 1.5 mm, H 6.0 mm PSA51	Shoulder Abutment 15° Ø 5.7 mm GH 3.0 mm, H 6.0 mm PSA52	Shoulder Abutment 15° Ø 5.7 mm GH 1.5 mm, H 4.0 mm PSA55	Shoulder Abutment 15° Ø 5.7 mm GH 3.0 mm, H 4.0 mm PSA56
Order No.				

				
Description	Shoulder Abutment 15° Ø 4.3 mm GH 1.5 mm, H 6.0 mm PSA61	Shoulder Abutment 15° Ø 4.3 mm GH 3.0 mm, H 6.0 mm PSA62	Shoulder Abutment 15° Ø 4.3 mm GH 1.5 mm, H 4.0 mm PSA65	Shoulder Abutment 15° Ø 4.3 mm GH 3.0 mm, H 4.0 mm PSA66
Order No.				


## Healing Abutments

				
Description	Healing Abutment Ø 6.4 mm GH 1.5 mm, Ø 5.7 mm	Healing Abutment Ø 6.4 mm GH 3.0 mm, Ø 5.7 mm	Healing Abutment Ø 5.0 mm GH 1.5 mm, Ø 4.3 mm	Healing Abutment Ø 5.0 mm GH 3.0 mm, Ø 4.3 mm
Order No.	PGF21	PGF22	PGF71	PGF72

## Occlusal Screw

	
Description	Occlusal Screw
Order No.	POS01




## Parallel Post




	
Description	Parallel Post for angled Abutments
Order No.	ZPS01




## TITANIUM BASE

### Abutments

			
Description	Titanium Base 0° GH 0.75mm, H 6.0 mm	Titanium Base 0° GH 1.5 mm, H 6.0 mm	Titanium Base 0° GH 3.0 mm, H 6.0 mm
Order No.	PTB01	PTB02	PTB03

			
Description	Titanium Base 0° GH 0.75 mm, H 4.0 mm	Titanium Base 0° GH 1.5 mm, H 4.0 mm	Titanium Base 0° GH 3.0 mm, H 4.0 mm
Order No.	PTB04	PTB05	PTB06

### Scan Base

	
Description	Scan Base
Order No.	PSB11

### Occlusal Screw



	
Description	Occlusal Screw - for Ceramics
Order No.	POSK1

#### Note:



The geometry of the occlusal screw - ceramic is designed such that shear stresses are avoided on the ceramic structure.




## MULTI UNIT ABUTMENT

### Abutments

		
Description	Multi Unit Abutment 0° GH 1.5 mm PMU01	Multi Unit Abutment 0° GH 3.0 mm PMU02
Order No.		

### Accessories




		
Description	Retention Coping short PRK01	Retention Coping long PRK02
Order No.		

			
Description	Occlusal Screw Multi Unit short POSM1	Occlusal Screw Multi Unit medium POSM2	Occlusal Screw Multi Unit long POSM3
Order No.			

			
Description	Gold Coping PSK01	Protective Cap PSKM1	Laboratory Analog one-piece PLAM1
Order No.			



	
Description	Screw Driver Multi Unit ZSDMU
Order No.	





Abutments




			
Description	LOCATOR® GH 2.0 mm PLOC2	LOCATOR® GH 3.0 mm PLOC3	LOCATOR® GH 4.0 mm PLOC4
Order No.			

Accessories

	
Description	LOCATOR Impression Copping
packed a	4
Order No.	PAKL1




				
Description	LOCATOR Female Analog Ø 4.0 mm	LOCATOR Male Processing Package	LOCATOR Male Processing Package	Black Processing Replacement Male
packed a	4	2	2	4
Order No.	PLAL1	PLLS1	PLLS2	PLPE1

				
Description	LOCATOR Block-Out Spacer	LOCATOR Processing Spacer	LOCATOR Abutment Retaining Sleeve	LOCATOR Parallel Post
packed a	20	4	4	4
Order No.	PLAR1	PLDH1	PLAH1	ZPP01





			
Description	LOCATOR Torque Wrench Driver	LOCATOR Torque Wrench Driver	LOCATOR Core Tool
Length	15 mm	21 mm	
Order No.	ZEIL1	ZEIL2	ZCT01

## BALL ANCHOR

### Abutments

			
Description	Ball Anchor GH 1.5 mm	Ball Anchor GH 3.0 mm	Ball Anchor GH 4.5 mm
Order No.	PKA01	PKA02	PKA03

### Accessories

				
Description	Laboratory Analog Ball Anchor	Matrix Ball Anchor 0.5 kg	Matrix Ball Anchor 1.0 kg	Matrix Ball Anchor 1.5 kg
Packed a	1	2	2	2
Order No.	PLAK1	PMAK1	PMAK2	PMAK3

				
Description	O-Ring red	O-Ring blue	O-Ring black	Seating Instrument Ball Anchor
Packed a	2	2	2	1
Order No.	PKOR1	PKOR2	PKOR3	ZEIK1


# SURGICAL INSTRUMENTS

## EXTENDED APICAL LENGTHS

When preparing the implant bed, one should take into account that the effective drilling depth varies apically by 0.4 - 0.6 mm and, depending on the diameter of the implant drill, is deeper than the desired implant length. This additional length must already be taken into account during the planning phase.


Tri-Spade Drill A	Tri-Spade Drill M	Tri-Spade Drill B
0.4 mm	0.5 mm	0.6 mm

### An example of extended apical length




## EXPLANATION OF LASER MARKINGS


### Twist drills




## Tri-Spade Drills




## Conical Reamer



## Cortical Countersink



## Tap



General Instruments




	Round Drill	Initial Bur short	Initial Bur long	Parallel Gauge
Description				
Order No.	ORB01	OIB01	OIB02	ORL01

Suitable for all  
implant Ø



	Twist Drill long	Twist Drill short	Drill Extension	High Torque Drill Extension
Description				
Order No.	OSB01	OSB02	ZBV01	ZBV02

Suitable for all  
implant Ø





## Tri-Spade Drill



A Ø 3.5 mm

Description	Tri-Spade Drill A XS	Tri-Spade Drill A S	Tri-Spade Drill A M
Order No.	ATS01	ATS02	ATS03



M Ø 4.0 mm

Description	Tri-Spade Drill M XS	Tri-Spade Drill M S	Tri-Spade Drill M M
Order No.	MTS01	MTS02	MTS03



B Ø 4.5 mm

Description	Tri-Spade Drill B XS	Tri-Spade Drill B S	Tri-Spade Drill B M
Order No.	BTS01	BTS02	BTS03

## Conical Reamer



Description

Conical Reamer  
A 6.6 mm / 8.0 mm

Conical Reamer

A 9.5 mm  
AKA02

Conical Reamer

A 11.0 mm  
AKA03

Conical Reamer

A 14.0 mm  
AKA04

Order No.

A Ø 3.5 mm



Description

Conical Reamer  
M 6.6 mm / 8.0 mm  
MKA01

Conical Reamer

M 9.5 mm  
MKA02

Conical Reamer

M 11.0 mm  
MKA03

Conical Reamer

M 14.0 mm  
MKA04

Order No.

M Ø 4.0 mm



Description

Conical Reamer  
B 6.6 mm / 8.0 mm  
BKA01

Conical Reamer

B 9.5 mm  
BKA02

Conical Reamer

B 11.0 mm  
BKA03

Conical Reamer




B 14.0 mm  
BKA04

Order No.




B Ø 4.5 mm


Note:  
This is an optional drill which can be used in cases of a very thick cortical plate to reduce mechanical pressure.

Cortical Countersink

			
Description	Cortical Countersink A	Cortical Countersink M	Cortical Countersink B
Order No.	AKS01	MKS01	BKS01




Tap

			
Description	Tap A (Ø 3.5 mm)	Tap M (Ø 4.0 mm)	Tap B (Ø 4.5 mm)
Order No.	AGS01	MGS04	BGS02



ACCESSORIES

			
Description	Ratchet	Torque Indicator	Torque Ratchet
Order No.	ZRA01	ZDMA1	ZRA02

			
Description	Open-end Wrench	High Torque Adapter S	High Torque Adapter M
Order No.	ZGS01	ZHT0S	ZHT0M

## Further Accessories



Description	Screw Driver S ZSD0S	Screw Driver L ZSD0L	Screw Driver short one-piece ZSD1S
Order No.			



Description	Drill Extension ZBV01	High Torque Drill Extension ZBV02
Order No.		



Description	Seating Instrument M ZEI0M	Seating Instrument L ZEI0L	Seating Instrument S Right Angle ZEI01	Seating Instrument L Right Angle ZEI02
Order No.				

### Note:

Check the O-ring before each use of the seating instruments. This must be replaced if damaged or worn. There is a risk that the grip between insertion instrument and implant will not be sufficient.



Description	Replacement O-Ring 10 ZEIOR
packed a Order No.	



Description	Abutment Remover	Abutment Remover S Right Angle	Abutment Remover L Right Angle
Order No.	ZAI01	ZAI51	ZAI52

**Note:**  
The self-locking cone enables an absolutely rotation-stable connection to the implant which provides a virtually bacteria-proof seal. This connection can be disengaged again with the aid of the abutment remover.



Description	Reverse Implant Driver
Order No.	ZRDI1



Description	Cover Screw Remover short	Cover Screw Remover long
Order No.	ZAD01	ZAD02




Description	Handwheel for Ratchet Insert	Handwheel	Handle for Laboratory
Order No.	ZHR01	ZHR02	ZHG01

## SURGICAL KITS

All instruments for surgical use are available in the myplant two surgical kit. This allows for particularly structured and user-friendly storage and the color coding of the instruments facilitates orientation during implant surgery.


- 1.** Torque Ratchet  
Order No. ZDMA1
- 2.** Ratchet  
Order No. ZRA01
- 3.** Open-end Wrench  
Order No. ZGS01
- 4.** Round Drill  
Order No. ORB01
- 5.** Initial Bur  
Order No. OIB01
- 6.** Initial Bur long  
Order No. OIB02
- 7.** Twist Drill  
Order No. OSB01
- 8.** Twist Drill short  
Order No. OSB02
- 9.** Parallel Gauge  
Order No. ORL01 (3x)



- 10.** Reverse Implant Driver  
Order No. ZRD11
- 11.** Handwheel  
Order No. ZHR01
- 12.** High Torque Adapter  
Order No. ZHT0S
- 13.** High Torque Adapter  
Order No. ZHT0M
- 14.** Drill Extension  
Order No. ZBV01
- 15.** High Torque Drill  
Extension  
Order No. ZBV02
- 16.** Screw Driver Right  
Angle short  
Order No. ZSD0S
- 17.** Screw Driver Right  
Angle long  
Order No. ZSD0L
- 18.** Screw Driver  
one-piece  
Order No. ZSD1S
- 19.** Seating Instrument  
Right Angle short  
Order No. ZE101



Description	Surgical Kit
Order No.	0CK01



Instruments for  
Implants with  
Ø 4.5 mm

Instruments for  
Implants with  
Ø 4.0 mm

Instruments for  
Implants with  
Ø 3.5 mm

Conical Reamer  
AØ 3.5: Order No. AKA01/ AKA02/ AKA03/ AKA04  
MØ 4.0: Order No. MKA01/ MKA02/ MKA03/ MKA04  
BØ 4.5: Order No. BKA01/ BKA02/ BKA03/ BKA04

Cortical Countersink  
Order No. AKS01/ MKS01/ BKS01

Tap  
Order No. AGS01 / MGS04 / BGS02

**20.** Seating Instrument  
Right Angle long  
Order No. ZEIO2  
**21.** Seating Instrument  
short  
Order No. ZEIO1M

**22.** Seating Instrument  
long  
Order No. ZEIO1L  
**23.** Cover Screw  
Remover short  
Order No. ZAD01

**24.** Cover Screw  
Remover  
Order No. ZAD02  
**25.** Abutment Remover  
Right Angle short  
Order No. ZAI51

**26.** Abutment Remover  
Right Angle long  
Order No. ZAI52  
**27.** Parallel Post  
Order No. ZPS01 (x2)


Surgical Kit (Order No. CKB40)  
contains the following articles:



Description Quantity / Kit Order No.	Seating Instrument M	Seating Instrument Right Angle	Screw Driver S	Reverse Implant Driver
	1	1	1	1
	ZEI0M	ZEI01	ZSD0S	ZRD11



Description Quantity / Kit Order No.	High Torque Adapter M	Drill Extension	Initial Bur short	Twist Drill	Open-end Wrench
	1	1	1	1	1
	ZHT0M	ZBV01	OIB01	OSB01	ZGS01



Description	Tri-Spade A XS	Tri-Spade A M	Tri-Spade M XS	Tri-Spade M M
Quantity / Kit	1	1	1	1
Order No.	ATS01	ATS03	MTS01	MTS03



Description	Conical Reamer M 6.6 mm / 8.0 mm	Conical Reamer M 9.5 mm	Conical Reamer M 11.0 mm	Conical Reamer M 14.0 mm	Tap M
Quantity / Kit	1	1	1	1	1
Order No.	MKA01	MKA02	MKA03	MKA04	MG504

# PROSTHETIC KITS

Prosthetic Kit (Order No. PPB01) contains the following articles:



Description	Torque Ratchet	Screw Driver one-piece	Screw Driver S	Screw Driver L	Cover Screw Remover short	Cover Screw Remover long
Quantity / Kit Order No.	1 ZRA02	1 ZSD1S	1 ZSD0S	1 ZSD0L	1 ZAD01	1 ZAD02




Description	High Torque Adapter M	Handwheel	Seating Instrument Ball Anchor	Positioning Key	Parallel Post
Quantity / Kit Order No.	1 ZHT0M	1 ZHR01	1 ZEIK1	1 ZAPS1	2 ZPS01

Laboratory Kit (Order No. PLB01)  
contains the following articles:




Description	High Torque Adapter M	Screw Driver S	Screw Driver L	Handwheel
Quantity / Kit	1	1	1	1
Order No.	ZHT0M	ZSD0S	ZSD0L	ZHR01




Drill-Stop-Control myplant two (Order No. BDSMP)  
contains the following articles:



Description	Twist Drill	Twist Drill	Twist Drill	Twist Drill
Ø in mm	020	024	029	033
Quantity / Kit	1	1	1	1
Order No.	TDS18	TDS18	TDS18	TDS18





Description	Drill Sleeve	Drill Sleeve	Drill Sleeve	Drill Sleeve
Length in mm	6.0	7.0	9.0	10.0
Quantity / Kit	1	1	1	1
Order No.	CL050	CL051	CL052	CL053

**Note:**  
All drill sleeves can be combined with the listed twist drills.



Description	Drill Sleeve	Drill Sleeve	Drill Sleeve	Drill Sleeve
Length in mm	11.0	12.0	13.0	14.0
Quantity / Kit	1	1	1	1
Order No.	CL054	CL055	CL056	CL057

## FURTHER INFORMATION


### INSTRUMENT PORTFOLIO

Taking the anatomical and spatial conditions into account, the suitable position and number of implants, as well as corresponding implant diameters and a suitable implant length are to be selected individually for every patient.

A systematic surgical technique for preparing the implant bed was developed for elevated primary stability (EPS). The implant bed can be prepared in three steps, adapted to the given bone quality. This leads to an optimization of the primary stability and expands the possibilities for immediate restorations.



















**Note:**

The instruments are supplied non-sterile. Prior to being used, the instruments must be checked for operational suitability. Please also observe the „General Instructions on Use and Safety Notes for MEISINGER products in the medical field“ and the „Notes on Reprocessing (cleaning, disinfection and sterilization) of medical devices by Hager & Meisinger GmbH“. Care should be taken to use the instruments intended specifically for the implant variant.
















## QUICK OVERVIEW OF SURGICAL PROCEDURE

	Ø 3.5 mm	Ø 4.0 mm	Ø 4.5 mm	opt. rotary speed	Torque
Smoothing of the Alveolar Ridge		 Round Drill 0RB01		2.000 min <sup>-1</sup>	-
Initial Drilling		 Initial Bur 0IB01 0IB02		1.000 min <sup>-1</sup>	-
Pilot Drilling Ø 2.0 mm		 Twist Drill 0SB01 0SB02		800 min <sup>-1</sup>	-
1. Expansion of the implant bed Ø 2.4 mm for A, M & B Implants		 Tri-Spade Drill A ATS01 ATS02 ATS03		800 min <sup>-1</sup>	-
2. Expansion of the implant bed Ø 2.9 mm for M & B Implants		 Tri-Spade Drill M MTS01 MTS02 MTS03		800 min <sup>-1</sup>	-
3. Expansion of the implant bed Ø 3.3 mm for B Implants			 Tri-Spade Drill B BTS01 BTS02 BTS03	800 min <sup>-1</sup>	-
Conical Expansion of the implant bed	 Conical Reamer A AKA01 / AKA02 / AKA03 / AKA04	 Conical Reamer M MKA01 / MKA02 / MKA03 / MKA04	 Conical Reamer B BKA01 / BKA02 / BKA03 / BKA04	50 min <sup>-1</sup>	max. 50 Ncm
Optional: Conical Extension of the Implant bed for low bone qualities and/or thick cortical plates	 Cortical Countersink A AKS01	 Cortical Countersink B MKS01	 Cortical Countersink C BKS01	50 min <sup>-1</sup>	max. 50 Ncm
Pretapping the Implant Thread	 Tap A AGS01	 Tap M MGS04	 Tap B BGS02	15 min <sup>-1</sup>	max. 50 Ncm
Implant Insertion	 Implant A A3566 / A3580 / A3595 / A3511 / A3514	 Implant M M4066 / M4080 / M4095 / M4011 / M4014	 Implant B B4566 / B4580 / B4595 / B4511 / B4514	15 min <sup>-1</sup>	max. 50 Ncm

## RECOMMENDED TORQUES FOR THE MYPLANT TWO PROSTHETICS

	Product	Torque
	Cover Screw	5 - 7 Ncm
	Healing Abutment	
	Transfer Post	
	Repositioning Post	
	Scan Base	
	Protective Cap	
	Occlusal Screw	10 Ncm
	Occlusal Screw for Ceramics	
	Occlusal Screw for Multi Unit	
	Standard Abutment	15 Ncm
	Shoulder Abutment	
	Titanium Base	
	Multi Unit Abutment	25 Ncm
	LOCATOR®	
	Ball Anchor	

**Note:** The images are examples only. The torque specifications are valid for all variants of the products mentioned.

## MATERIAL OVERVIEW

### UNALLOYED TITANIUM GRADE 4

**Titanium Grade 4 is used for implants and abutments.**

#### Chemical composition

O	Fe	C	N	H	Ti
0.4% max.	0.5% max.	0.08% max.	0.05% max.	0.015% max.	Rest

#### Technical notes

The titanium Grade 4 material complies with the standards ISO 5832-2 and ASTM F67.

### TITANIUM ALLOY GRADE 5-ELI / GRADE 23

**Titanium alloy Grade 5-ELI / Grade 23 is used for abutments.**

#### Chemical composition

Al	V	O	Fe	H	C	N	Ti
5.5-6.5% max.	3.5-4.5% max.	0.13% max.	0.25% max.	0.012% max.	0.08% max.	0.05% max.	Rest

#### Technical notes

The titanium alloy Grade 5-ELI material complies with the standards ISO 5832-3 and ASTM F136.

### STAINLESS STEEL 1.4197

**Stainless steel 1.4197 is used for instruments for implant bed preparation.**

#### Chemical composition

Cr	Mn	Mo	Si	Ni	C	S	P
12.5-14.0% max.	2.00% max.	1.00-1.50% max.	1.00% max.	0.75-1.50% max.	0.20-0.26% max.	0.15-0.27% max.	0.04% max.

#### Technical notes

The stainless steel 1.4197 material complies with the standards ISO 13504 and ASTM F899.

# IMPORTANT NOTES ON MYPLANT'S ORDERING SERVICES

## ORDER LINE AND SERVICE HOTLINE

We are at your service during the following times:

Monday – Thursday	8 a.m. – 6 p.m. Central European Time
Friday	8 a.m. – 5 p.m. Central European Time

## DELIVERY SERVICE

Deliveries within the European Union will take place within three days.

Information regarding deliveries outside the European Union can be provided upon request.

## RETURNS

- Only unused and unopened goods in original packaging.
- Within 30 days (valid from invoicing date) in conjunction with a copy invoice.

# MYPLANT TWO IMPLANT GUARANTEE

The myplant two guarantee from myplant GmbH includes a lifetime warranty on all implants in their product portfolio. This applies to manufacturing and material faults as well as implants which have not osseo-integrated. Implants will be replaced free of charge by Hager & Meisinger GmbH.

## EXCLUSION CRITERIA

Hager & Meisinger GmbH rejects any type of warranty or claims for damages due to inappropriate handling and non-compliance with the manufacturer's instructions.

Responsibility lies exclusively with the user and clinician.

## OTHER EXCLUSION CRITERIA

Loss of implant due to external causes (accidents, trauma, incorrect treatment, etc.).

Loss of implant due to non-compliance with contraindications (intake of medication, abuse of drugs and alcohol, consumption of tobacco, physical disorders). Detailed information is given in the instructions for use.

Additional claims and subsequent damages, for example, laboratory costs and clinical follow-up treatment are excluded from warranty.

## GUIDELINES

All myplant two items are to be used according to the instructions for use provided by the manufacturer. The use of non-original components as well as any type of modification may impair the function of the myplant two implant system and precludes any guarantee or replacement by Hager & Meisinger GmbH. This applies in particular to other application procedures that have not been recommended. System faults by mistaking tools and implants are to be minimized. Therefore, color coding and/or labelling should be observed. The processing and application of myplant two products takes place outside our control and is solely subject to the responsibility of the user. Any liability for such caused damages is excluded. Technical advice on the application of our products is verbal, in writing, via electronic media or demonstrations. This is based on state-of-the-art science and technology as known to us at the time of going to market. It does not absolve the user from the duty of personally checking the product, its indication-relevant suitability and training in the field of dento-alveolar surgery/implantology.

## AVAILABILITY

Some of the myplant two products given in this documentation may not be available in all countries. Detailed information is available upon request from myplant GmbH.

## ATTENTION

In addition to the warnings given in this document, our products are to be secured against aspiration when used intra-orally. Please observe the appropriate instructions for use as well as the manual for surgical and prosthetic procedures.

## VALIDITY


The publication of this document voids all previous versions.

## DOCUMENTATION

For further information on myplant two products, please refer to the relevant instructions for use or contact your myplant two representative.

## REGULATORY REQUIREMENTS

Meisinger stands for high quality medical devices since 1888. The quality management system of a company which manufactures medical devices must meet specific special requirements. These extremely high requirements are defined in ISO 13485 and meticulously complied with by our company. A MDSAP certificate according to ISO 13485:2016 confirms compliance with the requirements of international authorities in the USA (FDA), Canada (Health Canada), Australia (TGR), Japan (MHLW) and Brazil (ANVISA). All medical devices which you purchase from us as a customer, comply with all applicable requirements of the Medical Device Directive 93/42/EEC. Our company is certified by an independent Notified Body and certification is performed according to the specifications of standards. Current certificates can be found on our homepage [www.meisinger.de](http://www.meisinger.de)



## 1. Customer information

Name of physician: \_\_\_\_\_ Customer number: \_\_\_\_\_  
 Address: \_\_\_\_\_ Telephone: \_\_\_\_\_  
 \_\_\_\_\_ Email: \_\_\_\_\_  
 \_\_\_\_\_ Documented by: \_\_\_\_\_

## 2. Product information

REF No.:	LOT No.:	Date inserted:	Date removed:	Regio:
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

## 3. Patient information

Patient ID: \_\_\_\_\_ Age: \_\_\_\_\_ ☐ F ☐ M  
 Bone density ☐ D1 ☐ D2 ☐ D3 ☐ D4 Smoker? ☐ No ☐ Yes  
 Medical history:  
☐ Alcohol or drug abuse ☐ Blood coagulation disorder ☐ Chemotherapy during implantation  
☐ Diabetes mellitus ☐ Compromised immune resistance ☐ Treatment with corticosteroids  
☐ Lymphatic disorder ☐ Untreated endocrine disorders ☐ Psychic disorders  
☐ Radiation therapy in head/neck region ☐ Xerostomy ☐ No relevant findings  
☐ Immunological disorders ☐ Known allergies: \_\_\_\_\_  
☐ Other relevant disorders: \_\_\_\_\_

## 4. Surgical information

If the implant was inserted and removed on the same day, was another implant inserted successfully in the same place? ☐ No ☐ Yes LOT: \_\_\_\_\_  
 How was the implant inserted? ☐ Hand wheel ☐ Ratchet ☐ Angled handpiece Torque: \_\_\_\_\_ Ncm  
 Did problems occur with the pre-mounted transfer part? ☐ No ☐ Yes \_\_\_\_\_  
 Was one of the following points evident at the time of the intervention? ☐ Complication during preparation of the implant bed  
☐ Periodontal disease ☐ Mucosal disease ☐ Local infection / subacute chronic osteitis  
 What was the maximum speed employed during preparation? \_\_\_\_\_ min<sup>-1</sup>  
 Which drill was used last? \_\_\_\_\_  
 Was the thread tapped? ☐ Yes ☐ No  
 Was the enossal region covered completely by bone? ☐ Yes ☐ No  
 Was primary stability achieved? ☐ Yes ☐ No \_\_\_\_\_ Ncm  
 Was osseointegration achieved? ☐ Yes ☐ No  
 Was augmentation performed during surgery? ☐ Yes ☐ No Further information: \_\_\_\_\_  
 \_\_\_\_\_ Material used: \_\_\_\_\_  
 Was a membrane used? ☐ Yes ☐ No  
☐ Absorbable ☐ Not absorbable Material used: \_\_\_\_\_

## 5. Information about the event

Hygienic status around the implant? ☐ Very good ☐ Good ☐ Average ☐ Poor

Were one or more of the following factors involved in the event?

- |   |  |  |
|---|--|--|
| <input type="checkbox"/> Biomechanical overload | <input type="checkbox"/> Bruxism                               | <input type="checkbox"/> Bone resorption           |
| <input type="checkbox"/> Implant fracture       | <input type="checkbox"/> Peri-implantitis                      | <input type="checkbox"/> Infection                 |
| <input type="checkbox"/> Nerve compression      | <input type="checkbox"/> Sinus perforation                     | <input type="checkbox"/> Immediate implantation    |
| <input type="checkbox"/> Trauma or accident     | <input type="checkbox"/> Overheating of the bone               | <input type="checkbox"/> Insufficient bone quality |
| <input type="checkbox"/> Prior bone graft       | <input type="checkbox"/> Adjacent endodontically treated tooth |  |
| <input type="checkbox"/> Other: _____           |  |  |

The following was observed at implant loss

- |                                       |  |                                   |
|---------------------------------------|--|-----------------------------------|
| <input type="checkbox"/> Abscess      | <input type="checkbox"/> Asymptomatic          | <input type="checkbox"/> Bleeding |
| <input type="checkbox"/> Inflammation | <input type="checkbox"/> Increased sensitivity | <input type="checkbox"/> Fistula  |
| <input type="checkbox"/> Instability  | <input type="checkbox"/> Pain                  | <input type="checkbox"/> Swelling |
| <input type="checkbox"/> Numbness     | <input type="checkbox"/> Hypersensitivity      |                                   |

Had the implant already been prosthetically restored? ☐ No ☐ Yes (please answer Point 6)

What was the reason for implant loss in your opinion? \_\_\_\_\_

## 6. Information on the prosthetics

Type of restoration ☐ Crown ☐ Bridge ☐ Partial prosthesis, maxillary ☐ Partial prosthesis, mandibular  
☐ Full prosthesis, maxillary ☐ Full prosthesis, mandibular ☐ Other \_\_\_\_\_

When was the abutment placed?

| D | D | M | M | Y | Y |

Date of temporary restoration

| D | D | M | M | Y | Y |

Date of final restoration

| D | D | M | M | Y | Y |

Date of removal

| D | D | M | M | Y | Y |

Was a torque attachment used? ☐ Yes \_\_\_\_\_ Ncm ☐ No ☐ Not known

Were check-ups performed? ☐ Yes ☐ No

Case description: \_\_\_\_\_

## 7. Instruments

Approximate number of applications ☐ First time ☐ 2-5 ☐ 6-10 ☐ > 10

Method of cleaning ☐ Manual ☐ Ultrasonic ☐ Thermal disinfectant

Method of sterilization ☐ Autoclaving ☐ Dry heat ☐ Chemical autoclaving

Brief description of the event: \_\_\_\_\_

## 8. Confirmation

**All returned products are to be autoclaved and labelled as "sterile".**

Please add all the information necessary about the disputed products in this warranty form under consideration of the Hager & Meisinger GmbH warranty conditions and send this form including the autoclaved products and any X-rays back to Hager & Meisinger GmbH. Please use a padded bag for shipment - the loss of individual parts during shipment voids the warranty.

Date: \_\_\_\_\_

Signature of physician: \_\_\_\_\_



# INDEX

Order - No.		Description	Page	
IMPLANTS	IMPLANTS			
	A3511	Implant Ø 3.5 mm L 11.0 mm	12	
	A3514	Implant Ø 3.5 mm L 14.0 mm	12	
	A3566	Implant Ø 3.5 mm L 6.6 mm	12	
	A3580	Implant Ø 3.5 mm L 8.0 mm	12	
	A3595	Implant Ø 3.5 mm L 9.5 mm	12	
	M4011	Implant Ø 4.0 mm L 11.0 mm	12	
	M4014	Implant Ø 4.0 mm L 14.0 mm	12	
	M4066	Implant Ø 4.0 mm L 6.6 mm	12	
	M4080	Implant Ø 4.0 mm L 8.0 mm	12	
	M4095	Implant Ø 4.0 mm L 9.5 mm	12	
	B4511	Implant Ø 4.5 mm L 11.0 mm	12	
	B4514	Implant Ø 4.5 mm L 14.0 mm	12	
	B4566	Implant Ø 4.5 mm L 6.6 mm	12	
	B4580	Implant Ø 4.5 mm L 8.0 mm	12	
	B4595	Implant Ø 4.5 mm L 9.5 mm	12	
COVER SCREWS				
PVS00	Cover Screw 0 mm	13		
PVS01	Cover Screw 1 mm	13		
PVS02	Cover Screw 2 mm	13		
PROSTHETICS	HEALING ABUTMENT			
	PGF15	Healing Abutment Standard Ø 4.2 mm GH 1.5 mm	13	
	PGF30	Healing Abutment Standard Ø 4.2 mm GH 3.0 mm	13	
	PGF45	Healing Abutment Standard Ø 4.2 mm GH 4.5 mm	13	
	PGF21	Healing Abutment Shoulder Ø 6.4 mm GH 1.5 mm, Ø5.7 mm	20	
	PGF22	Healing Abutment Shoulder Ø 6.4 mm GH 3.0 mm, Ø5.7 mm	20	
	PGF71	Healing Abutment Shoulder Ø 5.0 mm GH 1.5 mm, Ø4.3 mm	20	
	PGF72	Healing Abutment Shoulder Ø 5.0 mm GH 3.0 mm, Ø4.3 mm	20	
	TRANSFER POSTS			
	PUP01	Transfer Post S open tray	14	
PUP02	Transfer Post L open tray	14		
REPOSITIONING POSTS				
PRP01	Repositioning Post, closed tray	14		
PRP02	Repositioning Post, narrow, closed tray	14		
IMPLANT ANALOG				
PLI01	Implant Analog	14		
TRANSFER CAPS				
PAK01	Transfer Cap H 4.0 mm	16		
PAK02	Transfer Cap H 6.0 mm	16		







PROSTHETICS



























Order - No. Description		Page	
LABORATORY ANALOGS			
PLA04	Laboratory Analog one-piece 0° H 4.0 mm	16	
PLA06	Laboratory Analog one-piece 0° H 6.0 mm	16	
PLA54	Laboratory Analog one-piece 15° H 4.0 mm	16	
PLA56	Laboratory Analog one-piece 15° H 6.0 mm	16	
HEALING CAPS			
PHK01	Healing Cap H 4.0 mm	16	
PHK02	Healing Cap H 6.0 mm	16	
STANDARD ABUTMENTS			
PAB01	Abutment 0° GH 1.5 mm, H 4.0 mm	15	
PAB02	Abutment 0° GH 1.5 mm, H 6.0 mm	15	
PAB03	Abutment 0° GH 3.0 mm, H 4.0 mm	15	
PAB04	Abutment 0° GH 3.0 mm, H 6.0 mm	15	
PAB05	Abutment 0° GH 4.5 mm, H 4.0 mm	15	
PAB06	Abutment 0° GH 4.5 mm, H 6.0 mm	15	
PAB51	Abutment 15° GH 1.5 mm, H 4.0 mm	15	
PAB52	Abutment 15° GH 1.5 mm, H 6.0 mm	15	
PAB53	Abutment 15° GH 3.0 mm, H 4.0 mm	15	
PAB54	Abutment 15° GH 3.0 mm, H 6.0 mm	15	
WAX-UP COPINGS			
PMK01	Wax-Up Coping H 4.0 mm	17	
PMK02	Wax-Up Coping H 6.0 mm	17	
PMK03	Wax-Up Coping H 4.0 mm Rotation Secured	17	
PMK04	Wax-Up Coping H 6.0 mm Rotation Secured	17	
TAPER CAPS			
PKK04	Taper Cap H 4.0 mm without Retention	17	
PKK06	Taper Cap H 6.0 mm without Retention	17	
PKK14	Taper Cap H 4.0 with Retention	17	
PKK16	Taper Cap H 6.0 with Retention	17	
TRIAL ABUTMENTS			
PAP00	Trial Abutment Kit	18	
PAP01	Trial Abutment 0° GH 1.5 mm, H 4.0 mm	18	
PAP02	Trial Abutment 0° GH 1.5 mm, H 6.0 mm	18	
PAP03	Trial Abutment 0° GH 3.0 mm, H 4.0 mm	18	
PAP04	Trial Abutment 0° GH 3.0 mm, H 6.0 mm	18	
PAP05	Trial Abutment 0° GH 4.5 mm, H 4.0 mm	18	
PAP06	Trial Abutment 0° GH 4.5 mm, H 6.0 mm	18	
PAP51	Trial Abutment 15° GH 1.5 mm, H 4.0 mm	18	
PAP52	Trial Abutment 15° GH 1.5 mm, H 6.0 mm	18	
PAP53	Trial Abutment 15° GH 3.0 mm, H 4.0 mm	18	
PAP54	Trial Abutment 15° GH 3.0 mm, H 6.0 mm	18	

Order - No.    Description			Page
SCHULTERABUTMENTS			
PSA01	Shoulder Abutment 0° Ø 5.7 mm GH 1.5 mm H 6.0 mm	19	
PSA02	Shoulder Abutment 0° Ø 5.7 mm GH 3.0 mm H 6.0 mm	19	
PSA05	Shoulder Abutment 0° Ø 5.7 mm GH 1.5 mm H 4.0 mm	19	
PSA06	Shoulder Abutment 0° Ø 5.7 mm GH 3.0 mm H 4.0 mm	19	
PSA11	Shoulder Abutment 0° Ø 4.3mm GH 1.5 mm H 6.0 mm	19	
PSA12	Shoulder Abutment 0° Ø 4.3mm GH 3.0 mm H 6.0 mm	19	
PSA15	Shoulder Abutment 0° Ø 4.3mm GH 1.5 mm H 4.0 mm	19	
PSA16	Shoulder Abutment 0° Ø 4.3mm GH 3.0 mm H 4.0 mm	19	
PSA51	Shoulder Abutment 15° Ø 5.7 mm GH 1.5 mm H 6.0 mm	19	
PSA52	Shoulder Abutment 15° Ø 5.7 mm GH 3.0 mm H 6.0 mm	19	
PSA55	Shoulder Abutment 15° Ø 5.7 mm GH 1.5 mm H 4.0 mm	19	
PSA56	Shoulder Abutment 15° Ø 5.7 mm GH 3.0 mm H 4.0 mm	19	
PSA61	Shoulder Abutment 15° Ø 4.3mm GH 1.5 mm H 6.0 mm	19	
PSA62	Shoulder Abutment 15° Ø 4.3mm GH 3.0 mm H 6.0 mm	19	
PSA65	Shoulder Abutment 15° Ø 4.3mm GH 1.5 mm H 4.0 mm	19	
PSA66	Shoulder Abutment 15° Ø 4.3mm GH 3.0 mm H 4.0 mm	19	
OCCLUSAL SCREWS			
POS01	Occlusal Screw	17 / 20	
POSK1	Occlusal Screw for Ceramics	21	
TITANIUM BASE AND EQUIPMENT			
PSB11	Scan Base	21	
PTB01	Titanium Base GH 0.75 mm, H 6.0 mm	21	
PTB02	Titanium Base GH 1.5 mm, H 6.0 mm	21	
PTB03	Titanium Base GH 3.0 mm, H 6.0 mm	21	
PTB04	Titanium Base GH 0.75 mm, H 4.0 mm	21	
PTB05	Titanium Base GH 1.5 mm, H 4.0 mm	21	
PTB06	Titanium Base GH 3.0 mm, H 4.0 mm	21	
MULTI UNIT ABUTMENTS AND EQUIPMENT			
PMU01	Multi Unit Abutment 0° GH 1.5 mm	22	
PMU02	Multi Unit Abutment 0° GH 3.0 mm	22	
PLAM1	Laboratory Analog Multi Unit	22	
POSM1	Occlusal Screw Multi Unit short	22	
POSM2	Occlusal Screw Multi Unit medium	22	
POSM3	Occlusal Screw Multi Unit long	22	
PRK01	Retention Coping Multi Unit short	22	
PRK02	Retention Coping Multi Unit long	22	
PSKM1	Protective Cap Multi Unit	22	
PSK01	Gold Coping Multi Unit	22	
ZSDMU	Screw Driver Multi Unit	22	


Order - No.    Description		Page		
PROSTHETICS	LOCATOR® AND ACCESSORIES			
	PLOC2	LOCATOR® Abutment GH 2.0 mm	23	
	PLOC3	LOCATOR® Abutment GH 3.0 mm	23	
	PLOC4	LOCATOR® Abutment GH 4.0 mm	23	
	PAKL1	LOCATOR® Impression Coping (4 Pcs)	23	
	PLAL1	LOCATOR® Female Analog (4 Pcs)	23	
	PLLS1	LOCATOR® Male Processing Package (2 Pcs)	23	
	PLLS2	LOCATOR® Male Processing Package - Extended Range (2 Pcs)	23	
	PLPE1	LOCATOR® Black Processing Replacement Male (4 Pcs)	23	
	PLAR1	LOCATOR® Block-Out Spacer (20 Pcs)	23	
	PLDH1	LOCATOR® Processing Spacer (4 Pcs)	23	
	PLAH1	LOCATOR® Abutment Retaining Sleeve (4 Pcs)	23	
	ZPP01	LOCATOR® Parallel Post (4 Pcs)	23	
	ZEIL1	LOCATOR® Torque Wrench Driver 15 mm	23	
	ZEIL2	LOCATOR® Torque Wrench Driver 21 mm	23	
	ZCT01	LOCATOR® Core Tool	23	
	BALL ANCHOR AND ACCESSORIES			
	PKA01	Ball Anchor GH 1.5 mm	24	
PKA02	Ball Anchor GH 3.0 mm	24		
PKA03	Ball Anchor GH 4.5 mm	24		
PLAK1	Laboratoty Analog Ball Anchor	24		
PMAK1	Matrix Ball Anchor 0.5 kg	24		
PMAK2	Matrix Ball Anchor 1.0 kg	24		
PMAK3	Matrix Ball Anchor 1.5 kg	24		
PKOR1	O-Ring red 0.5 kg	24		
PKOR2	O-Ring blue 1.0 kg	24		
PKOR3	O-Ring black 1.5 kg	24		
ZEIK1	Seating Insrument Ball Anchor	24 / 38		
INSTRUMENTS	DRILLS			
	ORB01	Round Drill	27	
	OIB01	Initial Bur short	27 / 36	
	OIB02	Initial Bur long	27	
	OSB01	Twist Drill long	27 / 36	
	OSB02	Twist Drill short	27	
	ZBV01	Drill Extension	27 / 32 / 36	
	ORL01	Parallel Gauge	27	
TRI-SPADE DRILLS				
ATS01	Tri-Spade Drill A XS	28 / 37		
ATS02	Tri-Spade Drill A S	28		
ATS03	Tri-Spade Drill A M	28 / 37		
MTS01	Tri-Spade Drill M XS	28 / 37		

INSTRUMENTS

Order - No.	Description	Page	
MTS02	Tri-Spade Drill M S	28	
MTS03	Tri-Spade Drill M M	28 / 37	
BTS01	Tri-Spade Drill B XS	28	
BTS02	Tri-Spade Drill B S	28	
BTS03	Tri-Spade Drill B M	28	
CONICAL REAMER			
AKA01	Conical Reamer A 6.6 mm / 8.0 mm	29	
AKA02	Conical Reamer A 9.5 mm	29	
AKA03	Conical Reamer A 11.0 mm	29	
AKA04	Conical Reamer A 14.0 mm	29	
MKA01	Conical Reamer M 6.6 mm / 8.0 mm	29 / 37	
MKA02	Conical Reamer M 9.5 mm	29 / 37	
MKA03	Conical Reamer M 11.0 mm	29 / 37	
MKA04	Conical Reamer M 14.0 mm	29 / 37	
BKA01	Conical Reamer B 6.6 mm / 8.0 mm	29	
BKA02	Conical Reamer B 9.5 mm	29	
BKA03	Conical Reamer B 11.0 mm	29	
BKA04	Conical Reamer B 14.0 mm	29	
CORTICAL COUNTERSINK			
AKS01	Cortical Countersink A	30	
MKS01	Cortical Countersink M	30	
BKS01	Cortical Countersink B	30	
TAPS			
AGS01	Tap A	30	
MGS04	Tap M	30 / 37	
BGS02	Tap B	30	
DRILL-STOP INSTRUMENTS			
TDS18	Twist Drill 2.0 mm Drill-Stop-Control	40	
TDS18	Twist Drill 2.4 mm Drill-Stop-Control	40	
TDS18	Twist Drill 2.9 mm Drill-Stop-Control	40	
TDS18	Twist Drill 3,3 mm Drill-Stop-Control	40	
CL050	Drill Sleeve 6.0 mm	41	
CL051	Drill Sleeve 7.0 mm	41	
CL052	Drill Sleeve 9.0 mm	41	
CL053	Drill Sleeve 10.0 mm	41	
CL054	Drill Sleeve 11.0 mm	41	
CL055	Drill Sleeve 12.0 mm	41	
CL056	Drill Sleeve 13.0 mm	41	
CL057	Drill Sleeve 14.0 mm	41	

Order - No.    Description			Page	
INSTRUMENTS	RATCHET AND RATCHET INSERTS			
	ZRA01	Ratchet	31	
	ZRA02	Torque Ratchet	31 / 38	
	ZDMA1	Torque Indicator	31	
	ZGS01	Open-end Wrench	31 / 36	
	ZHT0S	High Torque Adapter S	31	
	ZHT0M	High Torque Adapter M	31/36/38/39	
	ADDITIONAL ACCESSORY			
	ZSD0L	Screw Driver L	32/38/39	
	ZSD0S	Screw Driver S	32/36/38/39	
	ZSD1S	Screw Driver S (one-piece)	32 / 38	
	ZBV01	Drill Extension	27/32/36	
	ZBV02	High Torque Drill Extension	27 / 32	
	ZEI0M	Seating Instrument M	32 / 36	
	ZEI0L	Seating Instrument L	32	
	ZEI01	Seating Instrument short, Right Angle	32 / 36	
	ZEI02	Seating Instrument long, Right Angle	32	
	ZEIOR	Replacement O-Ring (10 Pcs)	32	
	ZHR01	Handwheel for Ratchet Insert	33/38/39	
	ZHR02	Handwheel	33	
	ZHG01	Handle Laboratory	33	
	ZAPS1	Positioning Key	38	
	ZPS01	Parallel Post	17/20/38	
	ZAI01	Abutment Remover	33	
	ZAI51	Abutment Remover Right Angle short	33	
	ZAI52	Abutment Remover Right Angle long	33	
	ZRD11	Reverse Implant Driver	33 / 36	
	ZAD01	Cover Screw Remover short	33 / 38	
	ZAD02	Cover Screw Remover long	33 / 38	
	ASSORTMENTS			
	0CK01	Surgical Kit	35	
	CKB40	Surgical Kit 4.0	36	
	BDSMP	Drill-Stop-Control myplant two	40	
	PPB01	Prosthetic Kit	38	
	PLB01	Laboratory Kit	39	





myplant GmbH | Hansemannstr. 10 | 41468 Neuss | Germany  
Phone: +49 2131 1259-465 | Fax: +49 2131 2012-222 | E-Mail: [info@myplant-dental.com](mailto:info@myplant-dental.com) | Internet: [www.myplant-dental.com](http://www.myplant-dental.com)



Hager & Meisinger GmbH | Hansemannstr. 10 | 41468 Neuss | Germany  
Phone: +49 2131 2012-0 | Fax: +49 2131 2012-222 | E-Mail: [info@meisinger.de](mailto:info@meisinger.de) | Internet: [www.meisinger.de](http://www.meisinger.de)

