

Screw System TX & Screw System TX Professional

Bone Fixation

Universal osteosynthesis screw set for easy and safe fixation of autologous bone cylinders, bone blocks and bone chips as part of pre-implantological augmentation.

NEW: With thread to screw head for bone augmentation via sandwich osteotomy!



Screw System TX (BTX00)




Screw System TX Professional (BTXPR)
Includes additional titanium screws in lengths
13.0 mm und 16.0 mm.



Our sales force would be pleased to
answer your questions: +49 2131 2012-303

Meisinger since
1888

Screw System TX

²Largest working part diameter in 1/10 mm

→ | ← External diameter ↗ ↘ Minimal diameter

Optionally available:

Art.-No. TL0T1

Screw System TX Professional

CE 0044

	Twist drills for TCT09		Twist drills for TCT10		Allport	Screw driver
Fig.	203RF	203RF	TC203	A2001	HM141F	31053
Shank ¹	205	205	204	204	205	205
Size ²	009	011	010	013	023	
Length in mm	9.0	9.0	12.0	12.0		

¹204=RA, 205=RA L

²Largest working part diameter in 1/10 mm



CE 0044

	Screws (7 Pcs. of each size)				Screws (3 Pcs. of each size)			
Fig.	TCT09*	TCT09*	TCT09*	TCT09*	TCT10*	TCT10*	TCT10*	TCT10*
Länge in mm	7.0	10.0	13.0	16.0	7.0	10.0	13.0	16.0
⌀*	0.9	0.9	0.9	0.9	1.0	1.0	1.0	1.0
□*	1.1	1.1	1.1	1.1	1.3	1.3	1.3	1.3

□* External diameter

⌀* Minimal diameter















With the reuse of disposable products, the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed

At a glance

- Universal screw set for easy and safe fixation of bone segments
- Safe retention and transmission of forces using insertion tool with internal hex connection for screws (Torx connection)
- Biocompatible screws if special titanium alloy with very flat screw head
- Very small screw diameter (0.90 mm) plus rescue screw (1.00 mm) for fixation of graft with high tensile strength
- Up to four different screw lengths (7, 10, 13, 16 mm)
- Optimal addition to all MEISINGER Bone Transfer-Systems and general augmentations

Fax for ordering +49 2131 2012-222

Screw System TX & Screw System TX Professional

	Figure	Lengths in mm	Minimal diameter	External diameter	Packaged à	Amount 
	BTX00	-	-	-	1	
	BTXPR	-	-	-	1	
	TCT09*	7.0	0.9	1.1	10	
	TCT09*	10.0	0.9	1.1	10	
	TCT09*	13.0	0.9	1.1	10	
	TCT09*	16.0	0.9	1.1	10	
	TCT10*	7.0	1.0	1.3	10	
	TCT10*	10.0	1.0	1.3	10	
	TCT10*	13.0	1.0	1.3	10	
	TCT10*	16.0	1.0	1.3	10	
	TL0T1	-	-	-	1	



CE 0044

Available from the dental trade.

Your trade depot:

Your delivery address:

Hager & Meisinger GmbH

Hansemannstr. 10 • 41468 Neuss • Germany
Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222
info@meisinger.de • www.meisinger.de

Meisinger USA, L.L.C.

10150 E. Easter Avenue • Centennial, Colorado 80112 • USA
Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
Toll free: +1 (866) 634-7464
info@meisingerusa.com • www.meisingerusa.com