

Prophy Twist Polisher

Our popular Twist Polishers: Now available for prophylaxis

The new Prophy Twist Polishers of MEISINGER convince through their effective and easy handling in the prophylaxis. The flexible polishing lamellae of the Prophy Twist Polishers adapt to any surface structure and generate flawless results on the natural tooth structure even in areas that are hard to reach, such as the occlusal surface. The disposable polishers have a plastic shank, thus ensuring a cost-effective application - a tangible value for every practice. The polishers are available in two sizes, each available as a test pack (10 pcs.) and a standard box (100 pcs.). For best results, we recommend using at a speed between 5,000 and 10,000 min⁻¹. The polishers can be used with and without polishing paste.



PTP01-140



PTP01-100

NEW

Disposable
instruments
for cost-effective
application!

At a glance

- Especially for polishing areas that are hard to reach
- Flawless results within a very short application time
- Flexible polishing lamellae adapt to any surface structure

Our sales force would be pleased to
answer your questions: +49 2131 2012-303

Meisinger
since
1888

Prophy Twist Polisher

Figure	Shank	Size	recomm. / max. revolutions	Packaged à
 PTP01 Prophy Twist Polisher	RA (204)	100	5.000 - 10.000 / 10.000 min ⁻¹	100
 PTP01 Prophy Twist Polisher	RA (204)	140	5.000 - 10.000 / 10.000 min ⁻¹	100
 PTP0T Prophy Twist Polisher (Testpackung)	RA (204)	100	5.000 - 10.000 / 10.000 min ⁻¹	10
 PTP0T Prophy Twist Polisher (Testpackung)	RA (204)	140	5.000 - 10.000 / 10.000 min ⁻¹	10

CE 0044 

Available from the dental trade.

Your trade depot:

Your delivery address:

Hager & Meisinger GmbH

Hansemannstr. 10 • 41468 Neuss • Germany
Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222
info@meisinger.de • www.meisinger.de

Meisinger USA. L.L.C.

10150 E. Easter Avenue • Centennial, Colorado 80112 • USA
Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
Toll free: +1 (866) 634-7464
info@meisingerusa.com • www.meisingerusa.com

85FL113E - 0118