

Crown Cutters

Deep X-cut for effective application

The new MEISINGER Crown Cutters have a stainless steel shank and a tungsten carbide working part. The working part is characterized by an extremely effective and deep X-cut. This toothing cuts absolutely reliable both on hart metals and on veneering ceramics, so it provides a broad field of application to the user. While using on veneering ceramics, working with these crown cutters is really time saving: A change of instruments isn't necessary while cutting at first the ceramic layer and then the inner part. For cutting of zirconia crowns we recommend the application of our Z838L.



HM4C



HM4CL



HM4CXL



HM4CXXL








Made in
Germany

At a glance

- Broad field of application
- Extremely effective, deep X-cut
- No change of instruments while working on veneering ceramics
- High service life
- Low wear

Our sales force would be pleased to
answer your questions: +49 2131 2012-303

Meisinger
since
1888

	Figur	Shank	Size	recommended / max. revolutions	Packaged à	Amount 
	HM4C	FG (314)	010	160.000 min ⁻¹ / 450.000 min ⁻¹	5	
	HM4C	FG (314)	012	160.000 min ⁻¹ / 450.000 min ⁻¹	5	
	HM4CL	FG (314)	010	160.000 min ⁻¹ / 450.000 min ⁻¹	5	
	HM4CL	FG (314)	012	160.000 min ⁻¹ / 450.000 min ⁻¹	5	
	HM4CXL	FG (314)	014	160.000 min ⁻¹ / 450.000 min ⁻¹	5	
	HM4CXXL	FG (314)	014	160.000 min ⁻¹ / 300.000 min ⁻¹	5	

CE 0044

Available from the dental trade.

Your trade depot:

Your delivery address:

Hager & Meisinger GmbH

Hansemannstr. 10 • 41468 Neuss • Germany
 Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222
 info@meisinger.de • www.meisinger.de

Meisinger USA. L.L.C.

10150 E. Easter Avenue • Centennial, Colorado 80112 • USA
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
 Toll free: +1 (866) 634-7464
 info@meisingerusa.com • www.meisingerusa.com