

# Occlu Twist Polishers

Precise polishing of the occlusal surface

The new Occlu Twist polishers from MEISINGER impress with their effective and intuitive application. The exceptional shape allows for efficient polishing of the occlusal surface and achieves flawless results even in deep fissures. This two-step polishing system produces high gloss results on all composite materials. Even hybrid composites, micro-filled composites and high-density restorative materials can be easily processed with the Occlu Twist polishers. Occlu Twist polishers are applicable without polishing paste and achieve brilliant polishing results in a very short processing time.



Made in  
Germany

## At a glance

- Specifically for occlusal surfaces
- Reaches deep fissures
- Intuitive application
- Brilliant polishing results
- Short processing time
- Extremely high service life

Our sales force would be pleased to  
answer your questions: +49 2131 2012-303

Meisinger  
since  
1888

# Occlu Twist Polishers

Fax for ordering +49 2131 2012-222

	Figure	Shank	Size	Recomm. / max. rotary speed	Packaged a	Amount 
	9765M	204	055	3.000 – 8.000 rpm / 20.000 rpm	2	
	9765F	204	055	3.000 – 8.000 rpm / 20.000 rpm	2	
For your further orders ...						

Available from the dental trade.

Your trade depot:

Your delivery address:

---



---



---



---



---



---



---



---

**Hager & Meisinger GmbH**

Hansemanstr. 10 • 41468 Neuss • Germany  
 Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222  
 info@meisinger.de • www.meisinger.de

**Meisinger USA. L.L.C.**

10150 E. Easter Avenue • Centennial, Colorado 80112 • USA  
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407  
 Toll free: +1 (866) 634-7464  
 info@meisingerusa.com • www.meisingerusa.com

85FL128EN - 04 19