

Ridge Contouring Bur

Surgical contouring and modification of the alveolar ridge

After tooth extraction or a longer period of edentulism, unwanted unevenness, rough surfaces or sharp edges may develop on the jawbone. In this case, the MEISINGER Ridge Contouring Bur offers the possibility of surgical contouring and modification of the alveolar ridge to smooth the bone. So, the alveolar ridge is prepared optimally for implant insertion or denture restoration.



At a glance



- Sharp and precise cutting edges for easy modification of the alveolar ridge
- Non toothed tip for optimal protection of the soft tissue
- High service life
- Reliable cutting performance

Our sales force would be pleased to
answer your questions: +49 2131 2012-303

Meisinger
since
1888

Ridge Contouring Bur

Fax for ordering +49 2131 2012-222

	Figure	Shank	Size	Max. rotary speed	Packaged a	Amount 
	RCB00	205	075	800 rpm	2	
For your further orders ...						

Available from the dental trade.

Your trade depot:	Your delivery address:
_____	_____
_____	_____
_____	_____
_____	_____

Hager & Meisinger GmbH
 Hansemannstr. 10 • 41468 Neuss • Germany
 Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222
 info@meisinger.de • www.meisinger.de

Meisinger USA, L.L.C.
 10150 E. Easter Avenue • Centennial, Colorado 80112 • USA
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
 Toll free: +1 (866) 634-7464
 info@meisingerusa.com • www.meisingerusa.com

85FL129EN - 0319