

# Ridge Contouring Bur

Surgical contouring and modification of the alveolar ridge

After tooth extraction or a longer period of edentulism, unwanted unevenness, rough surfaces or sharp edges may develop on the jawbone. In this case, the MEISINGER Ridge Contouring Bur offers the possibility of surgical contouring and modification of the alveolar ridge to smooth the bone. So, the alveolar ridge is prepared optimally for implant insertion or denture restoration.



## At a glance

- Sharp and precise cutting edges for easy modification of the alveolar ridge
- Non toothed tip for optimal protection of the soft tissue
- High service life
- Reliable cutting performance

Our sales force would be pleased to  
answer your questions: +49 2131 2012-303

*Meisinger* since 1888

