

CEREC® KIT

Digital Dental Academy

Universal system solution for your CAD/CAM restorations

We are able to present you with another MEISINGER highlight in the form of our new CEREC® Kit, developed in cooperation with international trainers at the Digital Dental Academy. This set lets you experience significant simplification and time-saving in your daily work. The attractive feature of the CEREC® Kit is its logical structure that guides you through your whole work procedure: Excavation – Preparation – Polish. Along with the most important instruments for the precise preparation of inlays, onlays, partial crowns and crowns, the kit also includes the outstanding MEISINGER crown separators HM7RX for metals and Z838L for zircon. The high-quality diamond polishers for a flawless finish complete the contents.





























At a glance

- Well-planned arrangement of instruments from preparation to high polish
- Precise preparation of inlays, onlays, partial crowns and crowns
- Perfect completion thanks to brilliant surface finish

Our sales force would be pleased to answer your questions: +49 2131 2012-303



CEREC® KIT DDA

| | Figure | Shank | Ref.-No. | Size | Packaged à | Amount  |
|---|--------|----------|-----------------|------|------------|--|
|  | DDA01 | - | - | - | - | |
|  | 808R | FG (314) | 806 314 238 524 | 012 | 5 | |
|  | HM7RX | FG (314) | 500 314 194 019 | 014 | 5 | |
|  | HM1SQ | RA (204) | 500 204 001 002 | 018 | 5 | |
|  | HM1SQ | RA (204) | 500 204 001 002 | 014 | 5 | |
|  | 850 | FG (314) | 806 314 198 524 | 021 | 5 | |
|  | 850F | FG (314) | 806 314 198 514 | 018 | 5 | |
|  | 850 | FG (314) | 806 314 198 524 | 016 | 5 | |
|  | 850F | FG (314) | 806 314 198 514 | 014 | 5 | |
|  | B878 | FG (314) | 809 314 298 544 | 018 | 5 | |
|  | B878 | FG (314) | 809 314 298 544 | 014 | 5 | |
|  | 830 | FG (314) | 806 314 257 524 | 023 | 5 | |
|  | 830F | FG (314) | 806 314 257 514 | 023 | 5 | |
|  | Z838L | FG (314) | 806 314 140 526 | 014 | 5 | |
|  | 863C | FG (314) | 806 314 250 504 | 012 | 5 | |
|  | 837L | FG (314) | 806 314 111 524 | 012 | 5 | |
|  | 837LF | FG (314) | 806 314 111 514 | 012 | 5 | |
|  | 859G | FG (314) | 806 314 166 534 | 014 | 5 | |
|  | 859F | FG (314) | 806 314 166 514 | 010 | 5 | |
|  | 9743M | RA (204) | 803 204 243 524 | 040 | 2 | |
|  | 9743F | RA (204) | 803 204 243 503 | 040 | 2 | |
|  | 9743C | RA (204) | 803 204 243 930 | 040 | 2 | |
|  | 9771M | RA (204) | 803 204 000 524 | 140 | 2 | |
|  | 9771F | RA (204) | 803 204 000 514 | 140 | 2 | |
|  | 9771C | RA (204) | 803 204 000 504 | 140 | 2 | |

CE 0044

Available from the dental trade.

Your trade depot:

Your delivery address:



Hager & Meisinger GmbH
 Hansemannstr. 10 • 41468 Neuss • Germany
 Phone: +49 2131 2012-0
 info@meisinger.de • www.meisinger.de

Meisinger USA. L.L.C.
 10200 E. Easter Avenue • Centennial, Colorado 80112 • USA
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
 Toll free: +1 (866) 634-7464
 info@meisingerusa.com • www.meisingerusa.com

Please observe the general instructions for the use and safety of MEISINGER products in the medical field as well as the notes on reprocessing (cleaning, disinfection and sterilization) Hager & Meisinger GmbH medical devices.

85FL122EN - 0817