

# PX Tungsten Carbide Cutters

## Effective processing of PEEK

The new MEISINGER PX Tungsten Carbide Cutters have a fine X-cut, which is distinguished by its optimized cutting angle. Due to that the new PX Tungsten Carbide Cutters are especially suitable for the treatment of difficult-to-machine PEEK. PEEK is a high-performance synthetic material that is used more and more frequently in dental implantology and prosthetics.

For processing it demands a high level of performance from the milling cutters in order to obtain clean results. The special toothing of the PX Tungsten Carbide Cutters allows for effective removal of this high-performance PEEK and care for perfect results free from streaks. In combination with a low contact pressure, the PX Tungsten Carbide Cutters provide smooth PEEK surfaces which are optimally prepared for final finishing.



HM73PX



HM138PX



HM129PX



HM79PX

Made in  
Germany





### At a glance

- Fine X-cut with an optimized cutting angle
- Effective processing of difficult-to-machine PEEK
- Processing of PEEK - free of streaks
- Provide smooth surfaces for final finishing
- Different working parts for optimal application on all restoration geometries
- High cutting performance

Our sales force would be pleased to  
answer your questions: +49 2131 2012-303

Meisinger  
since  
1888

# PX Tungsten Carbide Cutters Fax for ordering +49 2131 2012-222

	Figure	Shank	Size	recomm. / max. revolutions	Packaged à	Amount 
	HM73PX	HP (104)	014	15.000 min <sup>-1</sup> / 100.000 min <sup>-1</sup>	2	
	HM77PX	HP (104)	023	15.000 min <sup>-1</sup> / 100.000 min <sup>-1</sup>	2	
	HM79PX	HP (104)	040	15.000 min <sup>-1</sup> / 100.000 min <sup>-1</sup>	2	
	HM129PX	HP (104)	023	15.000 min <sup>-1</sup> / 100.000 min <sup>-1</sup>	2	
	HM138PX	HP (104)	023	15.000 min <sup>-1</sup> / 100.000 min <sup>-1</sup>	2	
	HM139PX	HP (104)	023	15.000 min <sup>-1</sup> / 100.000 min <sup>-1</sup>	2	
	HM261PX	HP (104)	023	15.000 min <sup>-1</sup> / 100.000 min <sup>-1</sup>	2	

CE 0044

Available from the dental trade.

Your trade depot:

Your delivery address:

---



---



---



---



---



---



---



---

**Hager & Meisinger GmbH**

Hansemannstr. 10 • 41468 Neuss • Germany  
 Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222  
 info@meisinger.de • www.meisinger.de

**Meisinger USA, L.L.C.**

10200 E. Easter Avenue • Centennial, Colorado 80112 • USA  
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407  
 Toll free: +1 (866) 634-7464  
 info@meisingerusa.com • www.meisingerusa.com

85FL117EN - 0317